

# Piccolo<sub>2</sub>-OEM

---

*FIBER-OPTIC OXYGEN METER  
OEM-MODULE WITH UART-INTERFACE*

*MANUAL*



The *Piccolo2-OEM* is manufactured by

***PyroScience*** GmbH  
Hubertusstr. 35  
52064 Aachen  
Germany

Phone	+49 (0)241 5183 2210
Fax	+49 (0)241 5183 2299
Email	<a href="mailto:info@pyroscience.com">info@pyroscience.com</a>
Internet	<a href="http://www.pyroscience.com">www.pyroscience.com</a>

Registered: Aachen HRB 17329, Germany

# TABLE OF CONTENT

<b>1</b>	<b>INTRODUCTION .....</b>	<b>1</b>
<b>2</b>	<b>SAFETY GUIDELINES .....</b>	<b>2</b>
<b>3</b>	<b>OVERVIEW .....</b>	<b>4</b>
3.1	OPTICAL PORT FOR OXYGEN SENSORS .....	4
3.2	EXTERNAL TEMPERATURE SENSOR .....	4
3.3	STATUS LED .....	6
<b>4</b>	<b>ELECTRICAL INTERFACE.....</b>	<b>7</b>
4.1	CONNECTOR .....	7
4.2	OPERATION VIA THE UART-INTERFACE .....	8
<b>5</b>	<b>SOFTWARE AND USB INTERFACE CABLE.....</b>	<b>9</b>
5.1	SOFTWARE INSTALLATION .....	9
5.2	USING THE PYRO OXYGEN LOGGER SOFTWARE.....	10
<b>6</b>	<b>TECHNICAL DRAWING.....</b>	<b>11</b>
<b>7</b>	<b>SPECIFICATIONS .....</b>	<b>12</b>

# 1 Introduction

The *Piccolo2-OEM* is a high precision optical oxygen meter for liquid or gas samples, which is based on the optoelectronics of the proven USB-device *Piccolo2* (item *PICO2*) from *PyroScience*. Most features of the *Piccolo2-OEM* are identical to the *Piccolo2*.

For a general introduction it is highly recommended to consult the detailed manual of the *Piccolo2*. Here we will cover mostly the differences of the *Piccolo2-OEM* compared to the *Piccolo2*.

The *Piccolo2-OEM* measures oxygen gas (e.g. %O<sub>2</sub>) or dissolved oxygen in liquids (DO). All sensors available for the *Piccolo2* can be used with the *Piccolo2-OEM*. The firmware provides internal calculation of all major oxygen units taking into account adjustable values of temperature, pressure, humidity, and salinity. In contrast to the *Piccolo2*, the *Piccolo2-OEM* includes additionally a temperature port, where a resistive temperature sensor (PT100, PT1000, or NTC; not included) can be directly connected. This temperature sensor can be used for automatic temperature compensation of the oxygen measurement. The communication interface is a serial interface with 3.0 V-levels (UART). An optional USB-interface-cable (item *PICO2-USB*) allows easy configuration, calibration and prototype testing of the module with the comfortable logger software *Pyro Oxygen Logger*. More information is available here:

Webpage of optical oxygen meter *Piccolo2*:

<https://www.pyroscience.com/piccolo2-optical-oxygen-meter.html>

Downloads for *Piccolo2* (manual, logger software):

<https://www.pyroscience.com/downloads.html>

## 2 Safety Guidelines

Before using the *Piccoloz-OEM* and its sensors, read carefully the instructions and user manuals.

In case of problems or damage, disconnect the instrument and mark it to prevent any further use! Consult *PyroScience* for advice! There are no serviceable parts inside the device. Please note that opening the housing will void the warranty!


The *Piccoloz-OEM* is not watertight. The *Piccoloz-OEM* should be kept under dry and clean conditions, avoiding moisture, dust, corrosive conditions and excessive heating of the instrument (e.g. direct sun light).

Calibration and application of the sensors is on the user's authority, as well as data acquisition, treatment and publication!

The sensors and the *Piccoloz-OEM* are not intended for medical or military purposes or any safety-critical applications.

The sensors should be used in the laboratory by qualified personnel only, following the user instructions and the safety guidelines of the manual, as well as the appropriate laws and guidelines for safety in the laboratory!

Keep the sensors and the *Piccoloz-OEM* out of reach of children!



**Figure 1: Overview Piccolo2-OEM**

## 3 Overview

Figure 1 provides an overview of the *Piccolo2-OEM*. The front provides the port for connecting fiber-optic oxygen sensors, as well as solder pads for an external temperature sensor (not included) enabling automatic temperature compensation of the oxygen measurement. The backside of the module provides the connector for the power supply and the digital communication interface, as well as a red status LED.

### 3.1 Optical Port for Oxygen Sensors


The *Piccolo2-OEM* provides a clamping screw for connecting optical fibers with an outer diameter of 3mm as known from the USB-device *Piccolo2*. Thus, all optical oxygen sensors available for the *Piccolo2* (e.g. robust probes, dipping probes, or “contactless” sensor solutions like sensor spots, flow through cells, respiration vials) can be directly connected to the *Piccolo2-OEM*. Please refer to the manual of the *Piccolo2* for further details.

### 3.2 External Temperature Sensor

The oxygen measurement is generally temperature dependent. Therefore, the *Piccolo2-OEM* offers automatic temperature compensation of the oxygen measurement, which typically should be activated if the temperature changes during measurements are more than ca. 1°C. The module offers for this a high-precision sensor interface which can be directly connected to a resistive temperature sensor (PT100, PT1000, or NTC; not included) which has to be placed within the sample where the oxygen is measured.

The firmware of the module offers special configuration registers which can be adapted to any PT100, PT1000 and NTC temperature sensor. A PT100- and PT1000-solution should be generally preferred, as they for most applications do not need any calibration, or in the worst case a simple offset calibration is

sufficient. An NTC sensor is generally more difficult to calibrate due to its non-linear temperature characteristics and its typ. large production tolerances.



**Figure 2: Connecting a resistive temperature sensor to the module**

The resistive temperature sensor has to be soldered to the 4 solder pads at the front of the module (Figure 2). For short distances (e.g. 10 cm) a simple 2-wire connection might be sufficient. For this it is important to shortcut the outer with the inner solder pads as indicated in Figure 2. For longer distances and/or for high precision measurements a 4-wire connection should be preferred (Figure 2).

**NOTE:** In order to minimize potential electrical noise coupling into the external temperature sensor, the cables should be twisted and as short as possible.



### 3.3 Status LED

The behavior of the status LED is given in Table 1.

Action/Status	Description	Behavior of status LED
Power-Up	The power supply is switched on.	A correct startup of the module is indicated by 4 flashes within 1-2 seconds, after this the LED remains lit.
Active	The module is either in idle mode waiting for a new command, or it is executing a command.	The LED is constantly lit.
Enter deep sleep	While the power supply is still enabled, the module can be put into deep sleep mode by the #STOP command.	The LED is switched off.
Exit deep sleep mode	The module reverts to active mode.	The LED is switched on.
#LOGO-command	The #LOGO-command is send to the communication interface	The LED flashes 4 times within 1-2 seconds, after this the LED remains lit.

**Table 1: Behavior of the status LED**

## 4 Electrical Interface

### 4.1 Connector

The electrical interface of the *Piccolo2-OEM* consists of the connector X1 (Figure 3). The package includes the fitting connector plug S1 (manufacturer: *Phoenix Contact*, type: *PTSMo,5/4-P-2,5*, item no.: *1778858*). Stripped cable ends can be connected to S1 without any soldering or crimping. When inserting or removing a stripped cable end (stripping length 6 mm, max. core diameter 0.5 mm<sup>2</sup>) into one of the connector holes of the connector S1, an internal spring mechanism has to be unlocked by pushing relatively strong with a small screw-driver (flat-bladed 2 mm in width) into the adjacent rectangular hole (Figure 3). The same manufacturer offers also fitting connector plugs for PCB mounting (details on request).




Figure 3: Electrical connectors of the Piccolo2-OEM

The pin configuration of the connector X1 is given in Table 2.

Pin	Name	Function	Description
1	VCC	Power	Power supply min. 3.3 VDC max. 5.0 VDC
2	RXD	Digital input 3.0 V levels (3.3 V & 5 V tolerant)	Data receive line of the UART interface
3	TXD	Digital output 3.0 V levels	Data transmission line of the UART interface
4	GND	Power	Ground

**Table 2: Pin configuration of the connector X1**

## 4.2 Operation via the UART-Interface

The *Piccolo2-OEM* is operated via a standard UART interface at 3.0 V levels (3.3 V and 5 V tolerant) consisting of a receive and a transmit line. The configuration of the UART-interface is as follows:

**19200 baud, 8 data bit, 1 stop bit, no parity, no handshake**

The baud rate can be changed by a special command of the communication protocol (refer to the communication protocol for further information). However, after power cycling the module always reverts to the standard baud rate of 19200.

Each action of the module is triggered by transmitting a text-command via the UART interface. For further details please refer to the communication protocol (available on request).

## 5 Software and USB interface cable

For the operation of the *Piccolo2-OEM* with a usual Windows PC, a coded USB interface cable (item *PICO2-USB*) is available from *PyroScience*. It includes a license for the comfortable oxygen logger software *Pyro Oxygen Logger*. Especially for initial testing purposes this logger software can speed up OEM-developments significantly. Additionally, the USB-interface cable *PICO2-USB* provides a virtual COM-port. Custom Windows-software can use this virtual COM-port for communicating directly with the module based on the communication protocol (available on request).

### 5.1 Software Installation

**System requirements:** PC with Windows 7/8/10 and min. 700 MB free disk space.

**IMPORTANT:** Do not connect the USB-interface cable *PICO2-USB* to your PC before the *Pyro Oxygen Logger* software has been installed. The software will install automatically the appropriate USB-drivers.

#### Installation steps:

- download the installer package and the manual for the newest version of the *Pyro Oxygen Logger* software from the *PyroScience* homepage:  
[www.pyroscience.com/downloads.html](http://www.pyroscience.com/downloads.html)
- unzip and start the installer and follow the instructions
- after the successful installation a new program group "Pyro Oxygen Logger" is added to the start menu, and a short-cut named "Oxygen Logger" can be found on the desktop
- connect one end of the USB-interface cable *PICO2-USB* to the connector X1 of the *Piccolo2-OEM*

- connect the other end to an USB-port of the PC. The status LED of the *Piccolo2-OEM* should flash shortly and remain lit, indicating the correct startup of the module.
- Start the *Pyro Oxygen Logger* software.

## 5.2 Using the Pyro Oxygen Logger software

The *Pyro Oxygen Logger* software with the USB interface cable *PICO<sub>2</sub>-USB* can be used for convenient configuration and calibration of the device. All device settings and calibration data are stored in non-volatile memory in the *Piccolo2-OEM* when the logger software is closed. Only minimal programming effort is necessary to trigger and read measurements from a device that is pre-configured this way. Thus, it is not necessary to implement any settings and calibrations routines in the custom software / firmware.

Please refer to the *Piccolo2*-manual for general operation instructions for the oxygen logger software.

## 6 Technical Drawing



The solder pads have 2.54mm pitch.

## 7 Specifications

General Specifications	
Dimensions	L=59 mm, Ø 17mm (without the optical port)
Weight	ca. 20 g
Power supply	min. 3.3 VDC max. 5.0 VDC
Power consumption	
-while acquiring a new oxygen data point	typ. 5-10 mA (typ. duration <100 ms)
-during idle time	typ. 2-3 mA
-during deep sleep mode	typ. 10 µA
Start-up time	
-from power off	1-2 s
-from deep sleep	ca. 200 ms
Interface	UART (3.0V levels, 5V tolerant), 19200 baud, 8 data bit, 1 stop bit, no parity, no handshake
Max. sample rate <sup>1</sup>	
-only with oxygen measurement	ca. 20 samples/s
-with oxygen and temperature measurement	ca. 5 samples/s
Operating temperature	0 to 50 °C
Storage temperature	-20 to 70 °C
Max. relative humidity	Non-condensing conditions

---

<sup>1</sup> Note: This max. sample rate refers only to the limits of the UART communication. It does not consider the actual response time of the connected optical oxygen sensor or of the temperature sensor.

<b>Oxygen Sensor</b>	
Oxygen measuring principle	lifetime detection of oxygen dependent REDFLASH indicator luminescence
Oxygen sensor connector	clamping screw for connecting optical fibers with an outer diameter of 3mm
Excitation wavelength	620 nm (orange-red)
Detection wavelength	760 nm (near infrared)
For further specifications of the oxygen sensors refer to the separately available specifications for the connected oxygen sensor.	
<b>Port for External Temperature Sensors</b>	
Compatible sensor types	PT100, PT1000, NTC (not included)
Measurement principle	2-wire or 4-wire resistance measurement via 24bit ADC
Resolution	<0.02 °C (typ. for PT100)
Accuracy	<+-0.2 °C (typ. PT100)
Range	-30 to 150 °C (typ. for PT100)
<b>Internal Temperature Sensor</b>	(located on internal PCB)
Resolution	0.02 °C
Accuracy	+ -0.3 °C
Range	-40 to 125 °C