


	<p>UN38.3 Lithium Battery Test Summary for PyroScience GmbH</p> <p>Products:</p> <ul style="list-style-type: none"> • FSGO2 • AquapHOx Logger
1	<p>Battery Manufacturer Spard New Energy Co., Ltd.</p>
2	<p>Product Manufacturer PyroScience GmbH Hubertusstrasse 35 52064 Aachen, Germany info@pyroscience.com +49 241 518 322 10</p>
3	<p>UN38.3 Test Lab Guangzhou MCM Certification & Testing Co., Ltd. 1 F No.13, Zhong San Section, Shi Guang Road Zhong Cun Street, Panyu District Guangzhou, Guangdong China</p>
4	<p>Test Report Number DYKJ20170724U02</p>
5	<p>Date of Test Report 25.09.2017</p>
6	<p>Description of Cell or Battery 3.7V, 1250mAh, 4.625Wh, 25g, lithium polymer pouch cell, Modell YT803448</p>
7	<p>UN38.3 Tests Performed and Successfully Passed Altitude simulation, thermal test, vibration, shock, external short-circuit, crush, overcharge, forced discharge</p>
8	<p>Assembled Battery Testing Requirements Not Applicable</p>
9	<p>Edition of UN Manual of Tests and Criteria Used UN Manual of Tests and Criteria ST/SG/AC.10/11/Rev.6, section 38.3</p>
10	<p>Name and Title of Signatory</p> <p></p> <p>Dr. Jan Fischer, PyroScience Engineer</p>