


Oxygen Sensor Integration Service

Customized Oxygen Sensor Minispots in your Device




- individual positioning of O₂ Sensor Minispots
- on glass, **COP/COC, PET, PMMA, PC, PS**
- spot diameter: 1 mm or 0.5 mm
- inside miniaturized structures

Integrate Oxygen Sensing into your Device!

Custom Integration of Oxygen Sensor Minispots

- biocompatible
- sterilizable (EtO, short-term 70% EtOH and IPA)
- compatible to individual processing steps (bonding methods, solvent and temperature exposure, etc.) - **catalogue available on request**
- **FirePlate-O₂** for High-Throughput Measurements

O₂



Make Oxygen visible in Your MPS!